

## Product Datasheet

### Norovirus Capsid Protein VP1 antibody (orb499859)

|                            |                                                                                             |
|----------------------------|---------------------------------------------------------------------------------------------|
| <b>Description</b>         | Norovirus Capsid Protein VP1 antibody                                                       |
| <b>Species/Host</b>        | Rabbit                                                                                      |
| <b>Conjugation</b>         | Unconjugated                                                                                |
| <b>Tested Applications</b> | ELISA                                                                                       |
| <b>Immunogen</b>           | KLH conjugated synthetic peptide derived from Norovirus Capsid protein VP1 (341-440/530 aa) |
| <b>Target</b>              | Norovirus Capsid protein VP1                                                                |
| <b>Preservatives</b>       | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.                       |
| <b>Form/Appearance</b>     | Liquid                                                                                      |
| <b>Concentration</b>       | 1mg/ml                                                                                      |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.            |
| <b>Note</b>                | For research use only                                                                       |
| <b>Isotype</b>             | IgG                                                                                         |
| <b>Clonality</b>           | Polyclonal                                                                                  |
| <b>MW</b>                  | 58 kDa                                                                                      |
| <b>Uniprot ID</b>          | <b>Q83884</b>                                                                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                             |