


Product Datasheet

Jak3 (phospho-Tyr785) antibody (orb499869)

Description	Jak3 (phospho-Tyr785) antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Jak3 around the phosphorylation site of Tyr785 (SD(p-Y)EL)
Target	Phospho-Jak3 (Tyr785)
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	125 kDa
Uniprot ID	P52333
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: Mouse Lung; Lane 2: Mo...