



Product Datasheet

Connexin-31 antibody (orb4999)

Description	Connexin-31 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Connexin-31 (51-150/270 aa)
Target	GJB3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	31 kDa
Uniprot ID	O75712
Expiration Date	12 months from date of receipt.



Western blot of rat myocardial tissue
us...



Immunohistochemical analysis of paraffin...