


Product Datasheet

UAP56 antibody (orb499907)

Description	UAP56 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human UAP56 (1-100/428 aa)
Target	UAP56
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	49 kDa
Uniprot ID	Q13838
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: Cerebrum (Mouse); Lane 2: Heart (Mouse); Lane 3: Cerebrum (Rat); Lane 4: Liver (Rat); Lane 5: AG-535 (Rat)