

Product Datasheet

Beta-Amyloid (1-40) antibody (orb499920)

Description	Beta-Amyloid (1-40) antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human beta-Amyloid (1-40/42 aa)
Target	beta-Amyloid (1-40)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	4.2 kDa
Uniprot ID	P05067
Expiration Date	12 months from date of receipt.