

Product Datasheet

Alpha-actinin-1 antibody (orb499943)

Description	Alpha-actinin-1 antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	Recombinant human Alpha-actinin-1 protein 1-205/894 (C-6x His-Tag)
Target	Alpha-actinin-1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	103 kDa
Uniprot ID	P12814
Expiration Date	12 months from date of receipt.