


Product Datasheet

Alpha Tubulin (Acetyl Lys40) antibody (orb499956)

Description	Alpha Tubulin (Acetyl Lys40) antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthesised acetylpeptide derived from human alpha Tubulin around the acetylation site of Acetyl Lys40 (SD(Ac-K)TI)
Target	alpha Tubulin (Acetyl Lys40)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	50 kDa
Uniprot ID	P68366
Dilution Range	WB=1:500-1000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: A549 (Human) Cell; La...