


Product Datasheet

Alpha Tubulin antibody (orb499957)

Description	Alpha Tubulin antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human alpha Tubulin Protein
Target	alpha Tubulin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	50 kDa
Uniprot ID	P68366
Dilution Range	WB=1:50000-200000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



Western blot:Lung
(Mouse);Hela(Human)
Ce...