

Product Datasheet

PPAR Delta antibody (orb499977)

Description	PPAR Delta antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-Fr, WB
Immunogen	Recombinant human PPAR delta Protein
Target	PPAR delta
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	48 kDa
Uniprot ID	P35396
Expiration Date	12 months from date of receipt.