


Product Datasheet

GFP-Tag antibody (orb499984)

Description	GFP-Tag antibody
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Green Fluorescence Protein
Target	GFP
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P42212
Dilution Range	WB=1:50000-500000
Expiration Date	12 months from date of receipt.



Transformed
(GFP or EV)
E. coli cells
we...