



Product Datasheet

AMACR antibody (orb500004)


Description	AMACR antibody
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human AMACR protein (full length)
Target	AMACR
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	42 kDa
Uniprot ID	Q9UHK6
Dilution Range	WB=1:500-1000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P analysis of Human kidney carcinoma...



Western blot:MCF-7(Human) Cell;HeLa(Huma...



Western blot:Lane 1:Mouse Kidney;Lane 2:...