


Product Datasheet

Human Activin A ELISA Kit (orb50001)

Description	Human Activin A ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: G311-S426
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and saliva.
MW	47442 MW
Uniprot ID	A4D1W7
Sensitivity	12 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 113-903 pg/ml.
Expiration Date	Please enquire.



Human
Activin A
ELISA Kit
standard
curve