

Product Datasheet

Beta II Tubulin antibody (orb500010)

Description	Beta II Tubulin antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human beta II Tubulin
Target	beta II Tubulin / Tubulin beta 2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	50 kDa
Uniprot ID	Q13885
Expiration Date	12 months from date of receipt.