


Product Datasheet

Mouse Activin A ELISA Kit (orb50002)

Description	Mouse Activin A ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: G311-S426
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and saliva.
MW	47392 MW
Uniprot ID	Q04998
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 62-422 pg/ml.
Expiration Date	Please enquire.



Mouse
Activin A
PicoKine
ELISA Kit
stand...