


Product Datasheet

Rat Activin A ELISA Kit (orb50003)

Description	Rat Activin A ELISA Kit
Reactivity	Rat
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: G311-S426
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and saliva.
MW	47406 MW
Uniprot ID	P18331
Sensitivity	12 pg/ml
Dilution Range	For serum, dilution ratio is 1:1 and concentration is 87-1759 pg/ml. For plasma, dilution ratio is 1:1 and concentration is 78-374 pg/ml.
Expiration Date	Please enquire.



Rat Activin A
PicoKine
ELISA Kit
standar...