


Product Datasheet

Human BDNF ELISA Kit (orb50004)

Description	Human BDNF ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: H129-R247
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA, citrate).
MW	27818 MW
Uniprot ID	P23560
Sensitivity	15 pg/ml
Dilution Range	For serum, dilution ratio is 1:100 and concentration is 7233-69396 pg/ml. For plasma, dilution ratio is 1:1 and concentration is 80-2694 pg/ml.
Expiration Date	Please enquire.



Line graph illustrates about the Ag-Ab r...