


Product Datasheet

Rat BDNF ELISA Kit (orb50005)

Description	Rat BDNF ELISA Kit
Reactivity	Rat
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: H129-R247
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA, citrate).
MW	28109 MW
Uniprot ID	P23363
Sensitivity	15 pg/ml
Dilution Range	Dilution ratio of 1:10, concentration in serum is around 5 ng/ml.
Expiration Date	Please enquire.



Rat BDNF
PicoKine
ELISA Kit
standard
cur...