


Product Datasheet

Mouse BDNF ELISA Kit (orb50006)

Description	Mouse BDNF ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: H129-R247
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA, citrate).
MW	28123 MW
Uniprot ID	P21237
Sensitivity	15 pg/ml
Dilution Range	Dilution ratio of 1:10, concentration in plasma is around 500 pg/ml.
Expiration Date	Please enquire.



Mouse BDNF
ELISA Kit
standard
curve