



Product Datasheet

STAT3 (phospho-Ser727) antibody (orb500068)


Description	STAT3 (phospho-Ser727) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human STAT3 around the phosphorylation site of Ser727 (PM(p-S)PR)
Target	phospho-STAT3 (Ser727)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	88 kDa
Uniprot ID	P40763
Dilution Range	WB=1:500-2000, IHC-P=1:100-200, ICC=1:50-100, Flow-Cyt=2ug/Test, IHC-F=1:100-200
Expiration Date	12 months from date of receipt.



Flow cytometry of A431.Primary Antibody ...



ICC analysis of A431 cell using phospho-...



IHC-P analysis of rat brain using phospho-...