


Product Datasheet

Human BMP-2 ELISA Kit (orb50007)

Description	Human BMP-2 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: Q283-R396
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	bone tissue, cell culture supernatants and serum.
MW	44702 MW
Uniprot ID	C8C060
Sensitivity	2 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum is around 50 pg/ml.
Expiration Date	Please enquire.



Human BMP-2 PicoKine ELISA Kit standard ...