

Product Datasheet

RAB3A antibody (orb500070)

Description	RAB3A antibody
Species/Host	Mouse
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant human RAB3A
Target	RAB3A
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	25 kDa
Uniprot ID	P20336
Dilution Range	WB=1:200-1000
Expiration Date	12 months from date of receipt.