



Product Datasheet

CK10 antibody (orb500076)


Description	CK10 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human Cytokeratin 10 protein, around 150-250aa
Target	CK10/Cytokeratin 10
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	59 kDa
Uniprot ID	P13645
Dilution Range	WB=1:500-2000, IHC-P=1:50-200
Expiration Date	12 months from date of receipt.



WB analysis of Skin (Mouse) Lysate at 40...



IHC-P image of mouse skin tissue using C...



IHC-P image of rat skin tissue using CK1...