



Product Datasheet

Collagen XVII antibody (orb500086)


Description	Collagen XVII antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	Recombinant human Collagen XVII protein (1300-1450/1497 aa)
Target	Collagen XVII
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/1ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	150 kDa
Uniprot ID	Q9UMD9
Dilution Range	WB=1:500-2000, IHC-P=1:50-200, ICC=1:50, IHC-F=1:50-200, IF=1:50-200
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
mouse skin
using Colla...



IHC-P
analysis of
human skin
using Colla...



IHC-P
analysis of
rat skin
using
Collage...