

Product Datasheet

Human CD40 Ligand/TNFSF5/CD40LG ELISA Kit (orb50011)

Description Human CD40 Ligand/TNFSF5/CD40LG ELISA Kit

Reactivity Human

Range 62.5 pg/ml - 4,000 pg/ml

Tested Applications ELISA

Immunogen Expression system for standard: E.coli; Immunogen sequence: E108-L261

Storage Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)

Note For research use only

Sample Types cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).


MW 29274 MW

Uniprot ID **P29965**

Sensitivity 15 pg/ml

Dilution Range For serum, dilution ratio is 1:5 and concentration is 0-13741 pg/ml. For plasma, dilution ratio is 1:1 and concentration is 0-166 pg/ml.

Expiration Date Please enquire.



Human
sCD40L
PicoKine
ELISA Kit
standard...