


Product Datasheet

Mouse EGF ELISA Kit (orb50013)

Description	Mouse EGF ELISA Kit
Reactivity	Mouse
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: N977-R1029
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and urine.
MW	133072 MW
Uniprot ID	P01132
Sensitivity	1 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 13 pg/ml.
Expiration Date	Please enquire.



Mouse EGF
PicoKine
ELISA Kit
standard
cu...