


Product Datasheet

Human Fibronectin ELISA Kit (orb50019)

Description	Human Fibronectin ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Target protein isolated from human plasma
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	262625 MW
Uniprot ID	P02751
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:10000, concentration in serum and plasma is around 15 ug/ml.
Expiration Date	Please enquire.



Human
Fibronectin
PicoKine
ELISA Kit
sta...