


Product Datasheet

Rat Fibronectin ELISA Kit (orb50020)

Description	Rat Fibronectin ELISA Kit
Reactivity	Rat
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Target protein isolated from rat plasma
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	272511 MW
Uniprot ID	P04937
Sensitivity	15 pg/ml
Expiration Date	Please enquire.



Rat
Fibronectin
PicoKine
ELISA Kit
stand...