


Product Datasheet

Human Fractalkine/CX3CL1 ELISA Kit (orb50021)

Description	Human Fractalkine/CX3CL1 ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q25-R339
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) , saliva, urine and human milk.
MW	42203 MW
Uniprot ID	A0N0N7
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum is 102-854 pg/ml.
Expiration Date	Please enquire.



Human
Fractalkine/CX3CL1
PicoKine ELISA ...