


Product Datasheet

Rat GDNF ELISA Kit (orb50024)

Description	Rat GDNF ELISA Kit
Reactivity	Rat
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S78-I211
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	23619 MW
Uniprot ID	A7UGJ1
Sensitivity	4 pg/ml
Expiration Date	Please enquire.



Rat GDNF
PicoKine
ELISA Kit
standard
cur...