

Product Datasheet

BCL-6 antibody (orb500274)

Description	BCL-6 antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	KLH conjugated synthetic peptide derived from human BCL-6 (676-704/706 aa)
Target	BCL-6
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/1ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	78 kDa
Uniprot ID	P41182
Expiration Date	12 months from date of receipt.