

Product Datasheet

RAB27B antibody (orb500281)

Description	RAB27B antibody
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant human RAB27B
Target	RAB27B
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	25 kDa
Uniprot ID	O00194
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.