

Product Datasheet

Bik (phospho-Thr33) antibody (orb500290)

Description	Bik (phospho-Thr33) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Bik around the phosphorylation site of Thr33 (GM(p-T)DS)
Target	phospho-Bik (Thr33)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	18 kDa
Uniprot ID	Q13323
Expiration Date	12 months from date of receipt.