


## Product Datasheet

# Human IFN Gamma/IFNG/Interferon Gamma ELISA Kit (orb50030)

|                            |                                                                                                                                                                                   |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Human IFN Gamma/IFNG/Interferon Gamma ELISA Kit                                                                                                                                   |
| <b>Reactivity</b>          | Human                                                                                                                                                                             |
| <b>Range</b>               | 15.6 pg/ml - 1,000 pg/ml                                                                                                                                                          |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                             |
| <b>Immunogen</b>           | Expression system for standard: E.coli; Immunogen sequence: Q24-Q166                                                                                                              |
| <b>Storage</b>             | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| <b>Note</b>                | For research use only                                                                                                                                                             |
| <b>Sample Types</b>        | cell culture supernatants and serum.                                                                                                                                              |
| <b>MW</b>                  | 19348 MW                                                                                                                                                                          |
| <b>Uniprot ID</b>          | <b>P01579</b>                                                                                                                                                                     |
| <b>Sensitivity</b>         | 2 pg/ml                                                                                                                                                                           |
| <b>Dilution Range</b>      | Dilution ratio of 1:1, concentration in serum is less than the lowest standard, 15.6 pg/ml.                                                                                       |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                   |



Human IFN  
gamma  
PicoKine  
ELISA Kit  
stand...