


Product Datasheet

Rat IGF-1 ELISA Kit (orb50033)

Description	Rat IGF-1 ELISA Kit
Reactivity	Rat
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: G49-A118
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	17079 MW
Uniprot ID	P08025
Sensitivity	5 pg/ml
Dilution Range	Dilution ratio of 1:500, concentration in serum and plasma is around 142-340 ng/ml.
Expiration Date	Please enquire.



Rat IGF-1
PicoKine
ELISA Kit
standard
cu...