



Product Datasheet

PTEN antibody (orb500350)

Description	PTEN antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human PTEN
Target	PTEN
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	44 kDa
Uniprot ID	P60484
Dilution Range	IHC-P=1:100-500, WB=1:500-2000
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
rat kidney
using PTEN ...



IHC-P
analysis of
mouse kidney
using PTE...