



Product Datasheet

HP1 Gamma antibody (orb500361)


Description	HP1 Gamma antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant human HP1 gamma Protein
Target	HP1 gamma
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	21 kDa
Uniprot ID	Q13185
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
human rectal
carcinoma...



IHC-P
analysis of
mouse liver
using HP1 ...



IHC-P
analysis of
rat brain
using HP1
ga...