

Product Datasheet

Histone H3 (mono Methyl K9) antibody (orb500365)

Description

Histone H3 (mono Methyl K9) antibody

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

KLH conjugated synthesised methylpeptide derived from human Histone H3 around the methylation site of mono methyl K9 (AR(mono Methyl-K)ST)

Target

Histone H3 (mono methyl K9)

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG1

Clonality

Monoclonal

MW

15 kDa

CAS Number

Histone H3.1

Uniprot ID


[P68431](https://www.uniprot.org/uniprot/P68431)

Dilution Range


WB=1:500-2000, IHC-P=1:100-500, IHC-F=1:100-500, ICC=1:100-500, IF=1:100-500, ELISA=1:1000-5000

Expiration Date


12 months from date of receipt.



IHC-P analysis of rat pancreas using His...



IHC-P analysis of mouse intestine tissue...



Western blot:U2os(Human) Cell;Primary:An...