


Product Datasheet

Mouse IGFBP-1 ELISA Kit (orb50037)

Description	Mouse IGFBP-1 ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A26-N272
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	29570 MW
Uniprot ID	P47876
Sensitivity	3 pg/ml
Dilution Range	Dilution ratio of 1:10, concentration in serum and plasma is 2737-4530 pg/ml.
Expiration Date	Please enquire.



Mouse
IGFBP-1
PicoKine
ELISA Kit
standar...