

Product Datasheet

MALAT1 antibody (orb500388)

Description	MALAT1 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide derived from human MALAT1 (1-45/57 aa)
Target	MALAT1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	6.3 kDa
Uniprot ID	Q9UHZ2
Expiration Date	12 months from date of receipt.