


Product Datasheet

Rat IL-1 Alpha/IL-1F1/IL1A ELISA Kit (orb50040)

Description	Rat IL-1 Alpha/IL-1F1/IL1A ELISA Kit
Reactivity	Rat
Range	4.69 pg/ml - 300 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: S115-S270
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and serum.
MW	30856 MW
Uniprot ID	P16598
Sensitivity	0.5 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum is less than the lowest standard, 12.5 pg/ml.
Expiration Date	Please enquire.



Rat IL-1
alpha
PicoKine
ELISA Kit
standa...