


Product Datasheet

Human IL-1 Beta/IL-1F2/IL1B ELISA Kit (orb50042)

Description	Human IL-1 Beta/IL-1F2/IL1B ELISA Kit
Reactivity	Human
Range	3.9 pg/ml - 250 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A117-S269
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	30748 MW
Uniprot ID	P01584
Sensitivity	0.15 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 12.5 pg/ml.
Expiration Date	Please enquire.



Human IL-1
beta
PicoKine
ELISA Kit
stand...