

Product Datasheet

Phospholipase D1 antibody (orb500423)

Description	Phospholipase D1 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Phospholipase D1 (601-700/1074 aa)
Target	Phospholipase D1/PLD1
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	124 kDa
Uniprot ID	Q13393
Expiration Date	12 months from date of receipt.