


Product Datasheet

Human IL-4/Interleukin-4 ELISA Kit (orb50048)

Description	Human IL-4/Interleukin-4 ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: H25-S153
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	17492 MW
Uniprot ID	P05112
Sensitivity	1.5 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 12.5 pg/ml.
Expiration Date	Please enquire.



Human IL-4
PicoKine
ELISA Kit
standard
C...