


Product Datasheet

Rat IL-4/Interleukin-4 ELISA Kit (orb50050)

Description	Rat IL-4/Interleukin-4 ELISA Kit
Reactivity	Rat
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: C25-S147
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and serum.
MW	16247 MW
Uniprot ID	P20096
Sensitivity	1 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum is less than the lowest standard, 15.6 pg/ml.
Expiration Date	Please enquire.



Rat IL-4
PicoKine
ELISA Kit
standard
cur...