

Product Datasheet

Phosphoserine antibody (orb500504)

Description	Phosphoserine antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Phosphoserine
Target	Phospho-Serine
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
CAS Number	407-41-0
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.