



Product Datasheet

Histone H3 (tri Methyl K4) antibody (orb500505)

Description	Histone H3 (tri Methyl K4) antibody
Species/Host	Mouse
Reactivity	Bovine, Drosophila, Human, Mouse, Porcine, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthesised methylpeptide derived from human Histone H3 around the methylation site of Tri Methyl K4 (RT(Tri Methyl-K)QT)
Target	Histone H3 (tri methyl K4)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	15 kDa
CAS Number	Histone H3.1
Uniprot ID	P68431
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
rat testis
using
Histo...



IHC-P
analysis of
Human
kidney
using His...