

Product Datasheet

Histone H3 (mono Methyl K79) antibody (orb500506)

Description

Histone H3 (mono Methyl K79)

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications
Immunogen

KLH conjugated synthesised methylpeptide derived from human Histone H3 around the methylation site of mono methyl K79 (DF(mono methyl K)TD)

Target

Histone H3 (mono methyl K79)

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG1

Clonality

Monoclonal

MW

15 kDa

CAS Number

Histone H3.1

Uniprot ID

P68431

Dilution Range

IHC-P=1:100-500

Expiration Date

12 months from date of receipt.