



Product Datasheet

GST Tag antibody (orb500513)

Description	GST Tag antibody
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant GST Protein
Target	GST tag
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Dilution Range	WB=1:30000-100000
Expiration Date	12 months from date of receipt.



Western blot:GST-Tagged fusion protein O...



Western blot:Recombinant GST Protein;Pri...