


Product Datasheet

Human IL-10/Interleukin-10 ELISA Kit (orb50055)

Description	Human IL-10/Interleukin-10 ELISA Kit
Reactivity	Human
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: S19-N178
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	20517 MW
Uniprot ID	P22301
Sensitivity	0.5 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 7.8 pg/ml.
Expiration Date	Please enquire.



Human IL-10 PicoKine ELISA Kit standard ...