



Product Datasheet

SMAD7 (phospho-Ser249) antibody (orb500554)


Description	SMAD7 (phospho-Ser249) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human SMAD7 around the phosphorylation site of Ser249 (SD(p-S)QL)
Target	Phospho-SMAD7 (Ser249)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	46 kDa
Uniprot ID	O15105
Dilution Range	WB=1:500-2000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



Paraformaldehyde-fixed, paraffin embedde...



Sample: 293T
(Human) Lysate at
30 ug Pri...



Paraformaldehyde-fixed, paraffin embedde...